



Human TERT DNA Copy Number Quantification qPCR Assay Kit (HTCQ)

Catalog #9018

100 reactions

Product Description

Telomerase reverse transcriptase (TERT) is a vital enzyme that maintains telomere length, crucial for cellular immortality and genome stability. The human TERT gene (hTERT) copy number often varies in cancer cells, serving as a key marker in oncogenesis studies. Increased hTERT copy number is linked to cancer development, tumor progression, poor prognosis, and genomic instability by enhancing telomerase activity. During cell immortalization, hTERT integration into the genome is common, making hTERT copy number measurement essential for distinguishing normal cells from immortalized ones.

ScienCell's Human TERT Copy Number Quantification qPCR Assay Kit (HTCQ) delivers a robust and accurate method to quantify hTERT copy number, driving insights into cancer development. The kit detects hTERT copy number in cell lines immortalized via hTERT overexpression and validates their consistency. The hTERT primer set amplifies a 100 bp-long region of the human TERT gene on chromosome 5. The single copy reference (SCR) primer set targets a 100 bp-long region on human chromosome 17. The reference DNA sample consists of a 1:1 ratio of hTERT copies to SCR copies and serves as a reference for calculating the ratio of hTERT to SCR of target samples.

Kit Components

| Cat # | Component | Quantity | Storage |
|-----------|---|----------|---------|
| MB6018a-1 | 2X GoldNStart TaqGreen qPCR master mix, 1 mL | 2 vials | -20°C |
| 9018a | hTERT primer set, lyophilized | 1 vial | -20°C |
| 9018b | Human single copy reference (SCR) primer set, lyophilized | 1 vial | -20°C |
| 9018c | Nuclease-free H ₂ O | 4 mL | 4°C |
| 9018d | Reference DNA sample (hTERT: SCR = 1:1) | 100 µL | -20°C |

Additional Materials Required (Materials Not Included in Kit)

| Component | Product Name |
|----------------------|---|
| DNA isolation kit | SpeedDNA Isolation Kit (ScienCell, Cat #MB6918) |
| genomic DNA template | Customers' samples |
| qPCR plate or tube | |

Quality Control

The hTERT and SCR primer sets are validated for specificity by qPCR with melt curve analysis and gel electrophoresis. Amplification efficiency is confirmed by template serial dilution (see Appendices 1 and 2).

Product Use

HTCQ is for research use only. It is not approved for human or animal use, or for application in clinical or *in vitro* diagnostic procedures.

Shipping and Storage

The product is shipped on dry ice. Upon receipt, store the GoldNStart TaqGreen qPCR master mix (Cat #MB6018a-1) in the dark at -20°C in a manual defrost freezer, the primers (Cat #9018a and 9018b) and the reference DNA sample (Cat #9018d) at -20°C in a manual defrost freezer, and the nuclease-free H₂O (Cat #9018c) at 4°C. Once thawed, do NOT refreeze GoldNStart TaqGreen qPCR master mix (Cat #MB6018a-1), and keep in the dark at 4°C or on ice at all times.

Procedures

Important: *Only use polymerases with hot-start capability to prevent possible primer-dimer formation. Only use nuclease-free reagents in PCR amplification.*

1. Prior to use, allow vials (Cat #9018a and #9018b) to warm to room temperature.
2. Centrifuge the vials at 1,500x g for 1 minute.
3. Add 200 µl nuclease-free H₂O (Cat #9018c) to hTERT primer set (lyophilized, Cat #9018a) to make hTERT primer stock solution. Aliquot as needed. Store at -20°C in a manual defrost freezer. Avoid repeated freeze-and-thaw cycles.
4. Add 200 µl nuclease-free H₂O (Cat #9018c) to SCR primer set (lyophilized, Cat #9018b) to make SCR primer stock solution. Aliquot as needed. Store at -20°C in a manual defrost freezer. Avoid repeated freeze-and-thaw cycles.
5. For the reference DNA sample (Cat #9018d), prepare two qPCR reactions, one with hTERT primer stock solution, and one with SCR primer stock solution. Prepare 20 µl qPCR reactions for one well as shown in Table 1.

Table 1.

| | |
|---|--------------|
| Reference DNA sample | 2 µl |
| Primer stock solution (hTERT or SCR) | 2 µl |
| 2X GoldNStart TaqGreen qPCR master mix (Cat #MB6018a-1) | 10 µl |
| Nuclease-free H ₂ O (Cat #9018c) | 6 µl |
| Total volume | 20 µl |

6. For each genomic DNA sample, prepare two qPCR reactions, one with hTERT primer stock solution, and one with SCR primer stock solution. Prepare 20 µl qPCR reactions for one well as shown in Table 2.

Table 2.

| | |
|---|--------------|
| Genomic DNA template (0.5 – 5 ng/μl) | 2 μl |
| Primer stock solution (hTERT or SCR) | 2 μl |
| 2X GoldNStart TaqGreen qPCR master mix (Cat #MB6018a-1) | 10 μl |
| Nuclease-free H ₂ O (Cat #9018c) | 6 μl |
| Total volume | 20 μl |

- Seal the qPCR reaction wells. Centrifuge the plates or tubes at 1,500x g for 15 seconds. For maximum reliability, replicates are highly recommended (minimum of 3).
- Refer to Table 3 for qPCR program setup. The 2X GoldNStart TaqGreen qPCR master mix (Cat #MB6018a-1) contains SYBR[®] Green as the reporter dye and does not contain a ROX passive reference dye. If the qPCR instrument being used has a "ROX passive reference dye" option, please deselect this option.

Note: The primary factors that determine optimal annealing temperature are the primer length and primer composition. Based on the properties of hTERT and SCR primer sets (Cat #9018a and #9018b), we highly recommend an annealing temperature of 62°C as shown in Table 3:

Table 3.

| Step | Temperature | Time | Number of cycles |
|----------------------|-------------------------------|------------|------------------|
| Initial denaturation | 95°C | 10 min | 1 |
| Denaturation | 95°C | 20 sec | 40 |
| Annealing | 62°C | 20 sec | |
| Extension | 72°C | 45 sec | |
| Data acquisition | Plate read | | |
| <i>Optional</i> | <i>Melting curve analysis</i> | | 1 |
| Hold | 20°C | Indefinite | 1 |

Quantification Method: Comparative $\Delta\Delta C_q$ (Quantification Cycle Value) Method

Note: Please refer to your qPCR instrument's data analysis software for data analysis. The method provided here serves as guidance for quick manual calculations.

- For hTERT, ΔC_q (hTERT) is the quantification cycle number difference of hTERT between the target and the reference DNA samples.

$$\Delta C_q (\text{hTERT}) = C_q (\text{hTERT, target sample}) - C_q (\text{hTERT, reference sample})$$

Note: the value of ΔC_q (hTERT) can be positive, 0, or negative.

- For single copy reference (SCR), ΔC_q (SCR) is the quantification cycle number difference of SCR between the target and the reference DNA samples.

$$\Delta C_q (\text{SCR}) = C_q (\text{SCR, target sample}) - C_q (\text{SCR, reference sample})$$

Note: the value of ΔC_q (SCR) can be positive, 0, or negative.

3. $\Delta\Delta Cq = \Delta Cq (\text{hTERT}) - \Delta Cq (\text{SCR})$
4. The ratio of hTERT to SCR of the target sample

$$= 2^{-\Delta\Delta Cq}$$
5. The hTERT copy number of the target sample

$$= 2 \times 2^{-\Delta\Delta Cq}$$

Example Calculations: Comparative $\Delta\Delta Cq$ (Quantification Cycle Value) Method

Table 3. Cq (Quantification Cycle) values of hTERT qPCR (hTERT) and single copy reference qPCR (SCR) obtained for the genomic DNA samples.

| <i>Primer set</i> | <i>Target sample</i> | <i>Reference sample</i> |
|-------------------|----------------------|-------------------------|
| hTERT | 24.83 | 20.99 |
| SCR | 25.64 | 20.80 |

$$\begin{aligned}
 \Delta Cq (\text{hTERT}) &= Cq (\text{hTERT, target sample}) - Cq (\text{hTERT, reference sample}) \\
 &= 24.83 - 20.99 \\
 &= 3.84
 \end{aligned}$$

$$\begin{aligned}
 \Delta Cq (\text{SCR}) &= Cq (\text{SCR, target sample}) - Cq (\text{SCR, reference sample}) \\
 &= 25.64 - 20.80 \\
 &= 4.84
 \end{aligned}$$

$$\begin{aligned}
 \Delta\Delta Cq &= \Delta Cq (\text{hTERT}) - \Delta Cq (\text{SCR}) \\
 &= 3.84 - 4.84 \\
 &= -1
 \end{aligned}$$


The ratio of hTERT to SCR of the target sample

$$\begin{aligned}
 &= 2^{-\Delta\Delta Cq} \\
 &= 2^{-(-1)} \\
 &= 2
 \end{aligned}$$


Conclusions: The average hTERT copy number in the target genomic DNA sample is twice that of the SCR control. Assuming there are 2 copies of SCR per diploid cell, this suggests there are 4 copies of hTERT per cell ($2 \times 2 = 4$).

Appendix 1: Quality assessment of hTERT primer set

A



B.



C.




Figure 1. Quality assessment of hTERT primer set. (A) qPCR amplification curves using serially diluted hTERT template. (B) Derivation of qPCR efficiency of hTERT primer set. (C) Separation of hTERT qPCR product by gel electrophoresis.




Appendix 2: Quality assessment of Single copy reference (SCR) primer set**A.****B.****C.**

Figure 2. Quality assessment of Single copy reference (SCR) primer set. (A) qPCR amplification curves using serially diluted SCR template. **(B)** Derivation of qPCR efficiency of SCR primer set. **(C)** Separation of SCR qPCR product by gel electrophoresis.