



Applied Biological Materials Inc

Tel: 1-866-757-2414
Email: info@abmGood.com
Website: www.abmGood.com

DNAfectin™2100

Store at 4°C

Cat. No.	Description	Quantity
G2100	DNAfectin™ 2100	1.0 mg/1.0 ml

Description

DNAfectin™2100 comprises of four unique formulations of polycations and liposomes, which will guarantee high transfection efficiency and low cytotoxicity for any cell type including primary cells.

Transfection Protocol

Use the following conditions as guidelines to transfect mammalian cells in a 6-well or 35mm dish format. For other culture vessels, please refer to Table 1.

1. **Adherent Cells:** 18 to 24 hours prior to transfection, seed cells at a density of $1-3 \times 10^5$ cells per well in 2.0 ml of appropriate growth medium (with serum and antibiotics if cells are cultured in the presence of them). Incubate the cells at 37°C in a CO₂ incubator until cells are 70% to 90% confluent at the time of transfection.

Suspension Cells: Just prior to preparing complexes, plate $3-5 \times 10^5$ cells in 0.8 ml of serum free medium without antibiotics.

Since transfection efficiency is sensitive to culture confluence, it is important to maintain a standard seeding protocol from experiment to experiment.

2. For each transfection sample, prepare the complexes as follows:

Solution A: Dilute 2.0 µg of DNA into 100 µl of serum-free, antibiotic-free medium.

Solution B: Vortex DNAfectin™2100 thoroughly prior use, then dilute 10-20 µl of DNAfectin™2100 in 100 µl serum-free, antibiotic-free medium.

Incubate Solution A and B at room temperature for 5 minutes.

3. Combine the solutions, mix gently to ensure uniform distribution and incubate for 20 minutes at room temperature. For suspension cells, go directly to step 5.

NOTE: Complexes are stable at room temperature for 3-5 hours.

4. **Adherent Cells ONLY:** Add 0.8ml of serum-free, antibiotic-free medium to DNAfectin™ 2100-DNA complex. Mix solution gently.

5. **Adherent Cells:** Remove growth medium from the cells and add 1.0ml of DNAfectin™ 2100-DNA solution to the each well containing cells.

Suspension Cells: Add 0.2 ml of the DNAfectin™ 2100-DNA solution into each well containing suspension cells in 0.8 ml serum-free, antibiotic-free medium.

- After 5-8 hours, remove transfection solution and add 2.0 ml of the appropriate growth medium (with serum and antibiotics) or add 0.1 ml of FBS directly into each vessel. Incubate the cells at 37°C in a CO₂ incubator for a total of 18-24 hours.
- To make stable cell lines: Passage cells at a 1:10 (or higher dilution) into fresh growth medium 24 hours post transfection. Selection medium can be added the following day if desired.

Optimizing Transfection for Specific Cell Lines

To achieve the maximum transfection efficiency and low cytotoxicity, optimize the transfection conditions by varying cell density along with DNA and DNAfectin™ concentrations. Optimal results have been observed when cells were 80-90% confluent and DNA(μg): DNAfectin™2100 (μl) ratios were 1:1 to 1:5.

Table 1: Reagent Quantities for Different Culture Vessels

Culture Vessel	Volume of plating medium per well	DNA(μg) in medium volume (μl)	DNAfectin™2100 in medium volume (μl)	Transfection medium vol.
24-well	500 μl	0.2-0.4 μg in 25 μl	2-4 μl in 25 μl	0.4 ml
12-well	1 ml	0.5-0.8 μg in 100 μl	5-8 μl in 100 μl	0.6 ml
6-well	2 ml	1.0-2.0 μg in 100 μl	10-20 μl in 100 μl	0.8 ml
35mm	2 ml	1.0-2.0 μg in 100 μl	10-20 μl in 100 μl	0.8 ml
60mm	5 ml	3.0-6.0μg in 500 μl	30-75 μl in 500 μl	2.4 ml
10-cm	10 ml	8.0-16.0 μg in 1.5ml	90-200 μl in 800μl	6.4 ml