

# Epithelial Cell Growth Supplement-2 (EpiCGS-2) Catalog Number: 4162

## **Product Description**

Epithelial Cell Growth Supplement-2 (EpiCGS-2) is a medium supplement designed for the optimal growth of human normal esophageal epithelial cells *in vitro*. It is a sterile, concentrated (100X) solution which contains growth factors, hormones, and proteins necessary for the culture of normal human epithelial cells. The supplement is formulated (quantitatively and qualitatively) to provide a defined and optimally balanced growth environment that maximally promotes the growth of normal human epithelial cells *in vitro*. The supplement is designed as an additive for epithelial cell medium-2 (EpiCM-2, Cat. No. 4111) and should be used in conjunction with that medium.

### Components

EpiCGS-2 is packaged in the quantity of supplement suited for a 500 ml bottle of epithelial cell medium-2. When a 500 ml bottle of epithelial cell medium-2 is supplemented with EpiCGS-2, the final concentrations of the supplement components per milliliter will be BSA 5  $\mu$ g, transferrin 5  $\mu$ g, BPE 50  $\mu$ g, insulin 2.5  $\mu$ g, FGF-2 5 ng, EGF 1ng, Epinephrine 500 ng, hydrocortisone 0.5  $\mu$ g, Prostaglandin E<sub>2</sub> 10<sup>-8</sup> M and T<sub>3</sub> 30 nM. <u>Note:</u> Due to BPE in the growth supplement, formation of lipoproteins can cause precipitates to be present; the color may vary with different lots.

### **Product Use**

<u>EpiCGS-2 is for research use only</u>. It is not approved for human or animal use, or for application in *in vitro* diagnostic procedures.

### Storage

Store the EpiCGS-2 at -20°C before adding to epithelial cell medium.

### Shipping

Dry ice.

### Prepare for use

Thaw <u>EpiCGS-2</u> at 37°C. Gently tilt the <u>EpiCGS-2</u> tube several times during thawing to help the contents dissolve. Make sure the contents of the supplement are completely dissolved into solution before adding to the medium. Rinse the bottle and tubes with 70% ethanol, and then wipe to remove excess. Remove the cap, being careful not to touch the interior threads with fingers. Add <u>EpiCGS-2</u> and other components (FBS and P/S solution) into basal medium in a sterile field, mix well and then the reconstituted medium is ready for use. Since several components of epithelial cell medium are light-labile, it is recommended that the medium not be exposed to light for lengthy periods of time. If the medium is stable for one month.

*Caution: If handled improperly, some components of the medium may present a health hazard. Take appropriate precautions when handling it, including the wearing of protective clothing and eyewear.* Dispose of properly.