



**Human Renal Epithelial Cell total RNA
(HREpiC total RNA)**

Catalog Number: 4125

RNA Preparation:	HREpiC total RNA was prepared from monolayer culture of Human Renal Epithelial Cells using Qiagen RNeasy kit. Total RNA concentration was determined by DU 800 UV/Visible spectrophotometer. An aliquot of 10 μ g was stored at -80°C .
Applications	RT-PCR, cDNA synthesis, RNA protection assay, northern blot analysis, purification of mRNA for library construction, PCR-select cDNA subtraction, and DNA micro array.
Storage Condition:	Store the RNA immediately at -80°C upon arrival. To ensure stability, we recommend aliquote the RNA into single-use fraction to avoid repeated freeze/thaw cycle.
Quality Control:	Purity: A_{260}/A_{280} ratio of 1.9-2.1 RNA integrity: intact 18S rRNA band and 28S rRNA band verified by a denaturing agarose gel, with a 2:1 ratio of 28S:18S. Residual DNA: not detected by RT-PCR.
Concentration	10 μ g (eighth 1 mg/ml or 0.5 mg/ml)
Shipping	Dry ice