



## Human Hepatocyte cDNA (HH cDNA)

Catalog Number: 5204

### Description

ScienCell human microglia cDNA is an excellent source of PCR-ready, first strand cDNA. It is ideal for use in a variety of gene discovery and analysis applications including:

- Amplification of genes of interest with specific or degenerate primes
- Analysis of cell-specific gene expression
- Characterization of alternately spliced mRNA

Ten micrograms of total RNA is primed with random primer oligonucleotides and reverse-transcribed using high-quality reagents and enzymes.

### Contents and Storage

ScienCell human microglia cDNA is supplied in a volume of 20  $\mu$ l for more than 20 reactions.

Store at  $-80^{\circ}\text{C}$ . ScienCell human microglia cDNA is stable for 6 months when stored properly.

### Product Use

HM cDNA is for research use only. It is not approved for human or animal use, or for application in *in vitro* diagnostic procedures.

### Shipping

Dry ice.