



**Human Dermal Fibroblast-neonatal Cell Pellet
(HDF-n Cell Pellet)**
Catalog #CP2310

Quantity: 5 million cells

Product Description

Human Dermal Fibroblast-neonatal Cell Pellet (HDF-n Cell Pellet) is prepared from early passage Human Dermal Fibroblasts-neonatal. Each pellet contains 5 million cells and can be used for a variety of applications including PCR, western blotting, genomic DNA library construction, and gene expression profiling. Cell pellets from ScienCell Research Laboratories are convenient and cost effective for researchers as it eliminates the need to acquire expensive tissues.

Applications

- Genomic DNA
- cDNA
- Total RNA
- microRNA
- Protein Lysates

Product Use

HDF-n Cell Pellet is for research use only. It is not approved for human or animal use, or for application in *in vitro* diagnostic procedures.

Storage

Store cell pellet immediately at -80°C upon arrival. Avoid freeze/thaw cycles.

Shipping

Dry ice.