



**Applied Biological Materials Inc.**

Telephone: 1-866-757-2414

Email: [info@abmgood.com](mailto:info@abmgood.com)

Website: [www.abmGood.com](http://www.abmGood.com)

---

## Cryopreservation Medium

<b>Cat.No.</b>	<b>Unit</b>
TM024	100ml

---

<b>Cat. No.</b>	TM024
<b>Name</b>	Cryopreservation Medium
<b>Description</b>	<b>abm's</b> Cryopreservation Medium is a ready-to-use complete formulation suitable for the cryopreservation of various types of mammalian cells. With its unique formula, this medium can be used to freeze cells cultured in PriGrow I, II, III and IV, as well as classical media such as DMEM, DMEM/F12, RPMI 1640 and Ham's F-12. This product contains 10% DMSO alongside additional optimized components that improve post-thaw cell recovery. Every lot of <b>abm's</b> Cryopreservation Medium is tested for pH, sterility and endotoxin.
<b>Unit</b>	100ml
<b>Sterility</b>	Pass
<b>Appearance</b>	Liquid
<b>QC</b>	Each lot has been fully tested and optimized for its ability to support in vitro cell culture. The testing parameters include: adherence rate (where applicable), typical morphology, and metabolism of the tested cells. Additionally, the media has been tested for sterility, pH, osmolality and endotoxin levels.
<b>Storage Condition</b>	Store at -20°C to 4°C and protect from light.
<b>Disclaimer</b>	<b>abm</b> products are intended for laboratory research purposes only, unless noted otherwise. They are not intended for use in humans. <b>abm</b> makes no warranties or representations as to the accuracy of the information on this site. Citations from literature and provided for informational purposes only. <b>abm</b> does not warrant that such information has been shown to be accurate.

---

**Caution:** This product is for research use only and is not intended for therapeutic or diagnostic applications.

Please contact a technical service representative for more information (1-866-757-2414).