



Applied Biological Materials Inc.

Telephone: 1-866-757-2414

Email: info@abmgood.com

Website: www.abmGood.com

Doxycycline

Cat.No.	Unit
G261	5 x 1.0 ml; 1.0 mg/ml

Cat. No.	G261
Name	Doxycycline
Unit	5 x 1.0 ml; 1.0 mg/ml
Category	Buffers & Chemicals
Description	Doxycycline is a member of the tetracycline antibiotics group, and is often used to treat bacterial infections. It is also an inhibitor of MMPs (matrix metalloproteinases). Doxycycline is utilized in tetracycline controlled transcriptional activation in order to regulate the expression of transgenes in cell cultures.
Concentration	1.0 mg/ml
Storage Condition	Store at -20°C.

Caution: This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information (1-866-757-2414).