



Applied Biological Materials Inc.

Telephone: 1-866-757-2414

Email: info@abmgood.com

Website: www.abmGood.com

PriCoat™ T75 Flasks, Filter Cap

Cat.No.	Unit
G299-T75	5 Flasks

Cat. No.	G299-T75
Name	PriCoat™ T75 Flasks, Filter Cap
Unit	5 Flasks
Category	Equipment and Lab Consumables
Description	<p>Features of PriCoat™ T75 Flasks, Filter Cap:</p> <p>Crystal clear virgin polystyrene</p> <p>Upper triangular and wide base shape provides better stability</p> <p>Volume graduations on both sides</p> <p>Vacuum plasma TC treatment for better cell attachment</p> <p>Vent cap with 0.22 µm hydrophobic filters for gas exchange and minimise risk of contamination</p> <p>Sterilized by E-Beam</p> <p>Non-pyrogenic</p>
Flask Volume	75 cm ²

Caution: *This product is for research use only and is not intended for therapeutic or diagnostic applications.*

Please contact a technical service representative for more information (1-866-757-2414).