



Applied Biological Materials Inc.

Telephone:1-866-757-2414

Email:info@abmgood.com

Website:www.abmGood.com

Penicillin/Streptomycin Solution

Cat.No.	Unit
G255	100ml

Cat. No.	G255
Name	Penicillin/Streptomycin Solution
Description	Penicillin- Streptomycin contains 10,000 units of penicillin (base), 10,000µg of streptomycin (base) per mL in WFI water. It is supplied as a 0.22 µm-filtered, 100X frozen liquid.
Unit	100ml
Usage	As a 100X concentration, it is designed to be supplemented as 1 ml per 99 ml of basal medium, yielding an iso-osmotic solution at final dilution.
Sterility	Pass
Appearance	Liquid
Storage Condition	This product is stable for 2 years at -20°C.

Caution: *This product is for research use only and is not intended for therapeutic or diagnostic applications.*

Please contact a technical service representative for more information (1-866-757-2414).