



Human Colonic Smooth Muscle Cell Genomic DNA
(HCSMC gDNA)
Catalog Number: 2949

Quantity: 5 µg

Product Description

ScienCell's genomic DNAs are isolated from early passage of over 100 different human primary cell types by Qiagen Allprep kit. The quality and purity of genomic DNA are tested by spectrophotometer and electrophoresis. A₂₆₀/A₂₈₀ is between 1.8 and 2.0. Our genomic DNAs are provided ready to use for SNP analysis, DNA methylation analysis, Southern Blotting, and PCR. Getting genomic DNAs from ScienCell will save you valuable research time to immediately identify genes of interest; no longer will you need to acquire tissues or cells and isolate genomic DNA by yourself.

Quality Control

- The A₂₆₀/A₂₈₀ ratio of the genomic DNA is measured to be 1.8 - 2.0 using Beckman Coulter DU 800 spectrophotometer.
- PCR detected alpha-actin gene.

Application

- PCR, SNP analysis and Southern Blotting
- Genomic DNA library construction
- Genomic DNA methylation study
- Profiling study in gene expression

Storage Conditions

Store genomic DNA at -80°C. Avoid repeated freeze/thaw cycles.