



**Applied Biological Materials Inc.**

Telephone: 1-866-757-2414

Email: [info@abmgood.com](mailto:info@abmgood.com)

Website: [www.abmGood.com](http://www.abmGood.com)

---

## pRetro-E1-HA

**Cat.No.**

E-516

**Unit**

10 µg

---

**Cat. No.**

E-516

**Name**

pRetro-E1-HA

**Unit**

10 µg


**Description**

Retroviral gene transfer is a widely used technique for efficient introduction of genes into the genome of any dividing cell type for stable expression. abm's Retro-Easy™ retroviral expression systems provide a quick and easy way of producing high titer recombinant retroviruses for efficient gene transduction. The Retro-Easy-HA™ Expression System is a complete retroviral expression system with the pRetro-E1-HA vector.

**Vector Map**

[pRetro-E1-HA](#) (click blue link to view)

**Vector Map**



**Storage Condition** -20°C or below.

**Note**

Not for resale without prior written consent of abm. This product is distributed for laboratory research only.

---

**Caution:** This product is for research use only and is not intended for therapeutic or diagnostic applications.

Please contact a technical service representative for more information (1-866-757-2414).