



Applied Biological Materials Inc.

Telephone: 1-866-757-2414

Email: info@abmgood.com

Website: www.abmGood.com

DTT

Cat.No.	Unit
G025	3 x 1.0 ml

Cat. No.	G025
Name	DTT
Unit	3 x 1.0 ml
Category	Buffers & General Chemicals
Description	DTT is a strong reducing agent that is usually used to reduce or "deprotect" thiolated DNA. It can also reduce the disulfide bonds in proteins and can block the formation of intramolecular and intermolecular bonds between the cysteine residues of protein. DTT is often mixed with a solution of DNA, and then removed by techniques such as filtration or chromatography.
Concentration	1M
Storage Condition	Store at 4°C

Caution: This product is for research use only and is not intended for therapeutic or diagnostic applications.

Please contact a technical service representative for more information (1-866-757-2414).