



Applied Biological Materials Inc.

Telephone: 1-866-757-2414

Email: info@abmgood.com

Website: www.abmGood.com

Polybrene

Cat.No.	Unit
G062	1.0 ml; 0.8 mg/ml

Cat. No.	G062
Name	Polybrene
Description	As a cationic polymer, polybrene increases the efficiency of infection in some cells with retrovirus (cell culture) by neutralizing the charge repulsion between virions and the cell surface. It can also be utilized in protein sequencing.
Unit	1.0 ml; 0.8 mg/ml
Sterility	Store at -20°C.
Storage Condition	Store at -20°C.

Caution: *This product is for research use only and is not intended for therapeutic or diagnostic applications.*

Please contact a technical service representative for more information (1-866-757-2414).