



Applied Biological Materials Inc.

Telephone: 1-866-757-2414

Email: info@abmgood.com

Website: www.abmGood.com

IPTG

Cat.No.	Unit
G116	1.0 g

Cat. No.	G116
Name	IPTG
Unit	1.0 g
Category	Reagents and Chemicals
Description	IPTG is a synthetic analogue of lactose, which works to inactivate the lac repressor and promote lactose utilization by inducing beta-galactosidase synthesis. IPTG is used to determine the lac phenotype in blue/white colony screening.
Storage Condition	Store at -20°C.

Caution: *This product is for research use only and is not intended for therapeutic or diagnostic applications.*

Please contact a technical service representative for more information (1-866-757-2414).