

Adeno-hTERT Antisense Virus

Cat.No.
G208

Unit
1.0 ml

Cat. No. G208

Name Adeno-hTERT Antisense Virus

Unit 1.0 ml


Description Standard Titer (10^6 pfu/ml) Recombinant Adenovirus expressing antisense hTERT gene.

Application Cell immortalization.

Expression System Type Adenovirus

Vector Map [pAdeno-hTERT-antisense](#)

Vector Map



Caution This product is distributed for laboratory research only.

Caution: This product is for research use only and is not intended for therapeutic or diagnostic applications.

Please contact a technical service representative for more information (1-866-757-2414).