


Retro-SV40T Virus

Cat.No.	Unit
G212	5 x 50 µl

Cat. No.	G212
Name	Retro-SV40T Virus
Unit	5 x 50 µl
Description	Recombinant Retrovirus (10^8 IU/ml) expressing the SV40 large and small T antigens.
Application	Cell immortalization.
Vector Map	pRetro-E2 SV40
Vector Map	



Caution This product is distributed for laboratory research only.

Caution: This product is for research use only and is not intended for therapeutic or diagnostic applications.

Please contact a technical service representative for more information (1-866-757-2414).