



Applied Biological Materials Inc.

Telephone: 1-866-757-2414

Email: info@abmgood.com

Website: www.abmGood.com

LB Agar

Cat.No.	Unit
G247	500 g

Cat. No.	G247
Name	LB Agar
Unit	500 g
Category	Reagents and Chemicals
Description	LB Agar is used for the cultivation of <i>E. coli</i> in molecular microbiology techniques. This medium provides nutrients for the bacteria so that it does not have to make its own, resulting in thriving colonies of <i>E. coli</i> . Additionally, LB Agar is useful since it does not preferentially grow one type of bacteria over another.
Storage Condition	Store under cool and dry conditions.

Caution: This product is for research use only and is not intended for therapeutic or diagnostic applications.

Please contact a technical service representative for more information (1-866-757-2414).