



**Applied Biological Materials Inc.**

Telephone:1-866-757-2414

Email:info@abmgood.com

Website:www.abmGood.com

---

## Gentamicin

<b>Cat.No.</b>	<b>Unit</b>
G272	10ml

---

<b>Cat. No.</b>	G272
<b>Name</b>	Gentamicin
<b>Description</b>	This product contains 50 mg/ml of active gentamicin sulfate in WFI water. It is 0.2 micron filtered and is tested for sterility.
<b>Unit</b>	10ml
<b>Usage</b>	Gentamicin has the advantage of being a broad-spectrum antibiotic that is non-toxic to animal cells and viruses at recommended concentrations. It should be added to cell cultures at a maximum concentration of 1ml per litre, with a working concentration of 50 µg/mL for eukaryotic cell culture and 15 µg/mL for prokaryotic cell culture.
<b>Sterility</b>	Pass
<b>Storage Condition</b>	This product is stored at 2-8°C. At room temperature, this product should be consumed within two days.

---

**Caution:** *This product is for research use only and is not intended for therapeutic or diagnostic applications.*

*Please contact a technical service representative for more information (1-866-757-2414).*