



Applied Biological Materials Inc.

Telephone: 1-866-757-2414

Email: info@abmgood.com

Website: www.abmGood.com

Antibiotic/Antimycotic (100X)

Cat.No.	Unit
G273	100ml

Cat. No.	G273
Name	Antibiotic/Antimycotic (100X)
Description	This product contains 10,000 units of penicillin G, 10,000µg of streptomycin sulfate, and 25µg of amphotericin B per ml in normal saline. It is 0.2 micron filtered and is tested for sterility.
Unit	100ml
Usage	This product is effective against gram positive bacteria, gram negative bacteria, and fungi. It should be added to cell cultures at a maximum concentration of 10ml per litre.
Sterility	Pass
Storage Condition	This product can be stored at -20°C for a maximum of 24 months. Avoid freeze-thaw cycles.

Caution: *This product is for research use only and is not intended for therapeutic or diagnostic applications.*

Please contact a technical service representative for more information (1-866-757-2414).