



**Applied Biological Materials Inc.**

Telephone:1-866-757-2414

Email:info@abmgood.com

Website:www.abmGood.com

---

## Amphotericin B (Fungizone)

<b>Cat.No.</b>	<b>Unit</b>
G274	100ml

---

<b>Cat. No.</b>	G274
<b>Name</b>	Amphotericin B (Fungizone)
<b>Description</b>	This product contains 250 µg/ml of amphotericin B and 205 µg/ml of deoxycholic acid in sterile ddH <sub>2</sub> O.
<b>Unit</b>	100ml
<b>Usage</b>	Amphotericin B is an antifungal agent used in cell cultures. It should be added to cell cultures at a maximum concentration of 10ml per litre of medium.
<b>Storage Condition</b>	This product is stored at -20°C. At 37°C, it is stable for up to 7 days. At 2-8°C, it is stable for up to 3 weeks. Avoid freeze-thaw cycles.

---

**Caution:** *This product is for research use only and is not intended for therapeutic or diagnostic applications.*

*Please contact a technical service representative for more information (1-866-757-2414).*