



Applied Biological Materials Inc.

Telephone:1-866-757-2414

Email:info@abmgood.com

Website:www.abmGood.com

PriCoat™ ECM T25 Flasks

Cat.No.	Unit
G999	10 Flasks

Cat. No.	G999
Name	PriCoat™ ECM T25 Flasks
Unit	10 Flasks
Category	Cell Culture Vessels
Description	Features of PriCoat™ ECM T25 Flasks: Coated with Extracellular Matrix to promote adhesion Crystal clear virgin polystyrene Upper triangular and wide base shape provides better stability Volume graduations on both sides Vacuum plasma TC treatment for better cell attachment Vent cap with 0.22 µm hydrophobic filters for gas exchange and minimise risk of contamination Sterilized by E-Beam Non-pyrogenic
Flask Volume	25 cm ²
Storage Condition	4°C to 30°C under dry conditions

Caution: *This product is for research use only and is not intended for therapeutic or diagnostic applications.*

Please contact a technical service representative for more information (1-866-757-2414).