



**Applied Biological Materials Inc.**

Telephone:1-866-757-2414

Email:info@abmgood.com

Website:www.abmGood.com

---

## Macrophage Culture Supplement/L-cell Supernatant

<b>Cat.No.</b>	<b>Unit</b>
TM142	40 ml

---

<b>Cat. No.</b>	TM142
<b>Name</b>	Macrophage Culture Supplement/L-cell Supernatant
<b>Description</b>	<b>abm's</b> Macrophage Culture Supplement/L-cell Supernatant may be used for cultivation of macrophages in culture. It is used for supplementation into media optimized for macrophage growth. Can be used in medium for up to 7 days.
<b>Unit</b>	40 ml
<b>Sterility</b>	Pass
<b>Appearance</b>	Liquid
<b>QC</b>	Each lot has been fully tested and optimized for its ability to support in vitro cell culture. The testing parameters include: adherence rate (where applicable), typical morphology, and metabolism of the tested cells. Additionally, the media has been tested for sterility, pH, osmolality and endotoxin levels.
<b>Storage Condition</b>	Store at -20°C.
<b>Disclaimer</b>	<b>abm</b> products are intended for laboratory research purposes only, unless noted otherwise. They are not intended for use in humans. <b>abm</b> makes no warranties or representations as to the accuracy of the information on this site. Citations from literature and provided for informational purposes only. abm does not warrant that such information has been shown to be accurate.

---

**Caution:** This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information (1-866-757-2414).